

Calibration Standards

CALIBRATION LABORATORY

Glassco's German calibration machines enables 0.003% accuracy which provides an A-class tolerance. Our automatic quality control center provides verification of Burettes, Pipettes & Flasks. All the Volumetric Glassware comes with Different kind of calibration Certificates.



QUALITY CONTROL **CENTRE MACHINE**



LARGE VOLUME CALIBRATION LAB MACHINE



SMALL VOLUMETRIC CALIBRATION LAB MACHINE



PIPETTE CALIBRATION MACHINE



National Accreditation Board for Testing and Calibration Laboratories (A Constituent Board of Quality Council of India)



CERTIFICATE OF ACCREDITATION

CALIBRATION LABORATORY, GLASSCO LABORATORY EQUIPMENTS PVT. LTD. has been assessed and accredited in accordance with the standard

ISO/IEC 17025:2005

"General Requirements for the Competence of Testing & Calibration Laboratories"

for its facilities at

P.O. Khudda Kalan, Manglai, Ambala Cantt., Haryana

in the field of

CALIBRATION



continued satisfactory compliance to the above standard & the relevant requirements of NABL through out ANAR. we (To see the maps of accordinator of this information, someone after

Signed for and on behalf of NAIL









Chief Exercitive Officer

Glassco's Calibration Laboratory has undertaken accreditation for its system for calibration and testing as

per ISO 17025:2005. Glassco's Calibration Laboratory is NABL approved. NABL has a MRA under the Asia Pacific Laboratory Accreditation Cooperation (APLAC) and the International Laboratory Accreditation Cooperation (ILAC), which means the calibration provided to customers for volumetric ware is acceptable in

calibration laboratories worldwide. Results of calibration are guaranteed for precision as per the DIN, ISO, USP & ASTM standards.

The calibrations performed by Glassco calibration laboratory are essential done under ISO 17025:2005.

Please visit our website www.glasscolabs.com to see list of accreditation bodies.