

PRODUCT DATA SHEET



CV0006

Screw Neck Vial, Squat Form, With Unattached PP Closure



Product comes with unattached PP Screw Cap

Specifications

Vial Material	Glass Borosilicate 3 (Autoclavable)
Capacity	14ml
Outer Diameter	23mm
Height	58mm
Closure Type	Screw, Black Polypropylene
Closure Liner	Aluminium
Pack Quantity	100
Packaging	Cardboard Box & Plastic Bag
Pack Length	385mm
Pack Width	220mm
Pack Height	67mm
Pack Weight	1.5Kg
Case Quantity	4packs (400vials)
Case Length	400mm
Case Width	240mm
Case Height	290mm
Case Weight	6.1Kg
Pallet Quantity	200packs (20000vials)

What's in the box?



Streamline your workflow and simplify your storage with our pack of 100 vials and caps. They neatly stack and nest together, maximizing your space and allowing convenient access whenever needed.

Unrivalled Quality

Each vial is crafted using premium-quality, chemically-resistant glass that ensures excellent durability and protection, and providing unparalleled clarity.

Range of Sizes

We understand that every project is unique. That's why our new range is available in a variety of sizes, from 2 ml to 100 ml, catering to all your specific storage needs.

Environmental Responsibility

We're committed to protecting the planet. Our unique and innovative packaging is 100% recyclable and eco-friendly, without compromising on quality or durability.

If you require any more information on this product, please contact our customer service team.