

PRODUCT DATA SHEET



CV0011

Specimen Tubes, With Unattached Push In, Polythene Stoppers



Product comes with unattached Push In, Polythene Stoppers

Specifications

Vial Material	Glass, Soda <i>(Not recommended for autoclave)</i>	●
Outer Diameter	10mm	
Height	50mm	
Closure Type	Push In, Polythene Stoppers	
Closure Liner	Unlined	
Pack Quantity	100	
Packaging	Cardboard Box & Plastic Bag	

What's in the box?



Streamline your workflow and simplify your storage with our pack of 100 vials and caps. They neatly stack and nest together, maximizing your space and allowing convenient access whenever needed.

Unrivalled Quality

Each vial is crafted using premium-quality, chemically-resistant glass that ensures excellent durability and protection, and providing unparalleled clarity.

Range of Sizes

We understand that every project is unique. That's why our new range is available in a variety of sizes, from 10mm to 100mm, catering to all your specific storage needs.

Environmental Responsibility

We're committed to protecting the planet. Our unique and innovative packaging is 100% recyclable and eco-friendly, without compromising on quality or durability.

If you require any more information on this product, please contact our customer service team.