

PRODUCT DATA SHEET



CV0013

Screw Neck Vial, Tall Form, With Unattached PP Closure



Product comes with unattached PP Screw Cap

Specifications

Vial Material	Glass Borosilicate 3 (Autoclavable)
Capacity	1.75ml
Outer Diameter	11mm
Height	36mm
Closure Type	Screw, Black Polypropylene
Closure Liner	Aluminium
Pack Quantity	100
Packaging	Cardboard Box & Plastic Bag
Pack Length	205mm
Pack Width	140mm
Pack Height	40mm
Pack Weight	0.4Kg
Case Quantity	16packs (1600vials)
Case Length	430mm
Case Width	310mm
Case Height	190mm
Case Weight	6.6Kg
Pallet Quantity	1000packs (100000Vials)

What's in the box?



Streamline your workflow and simplify your storage with our pack of 100 vials and caps. They neatly stack and nest together, maximizing your space and allowing convenient access whenever needed.

Unrivalled Quality

Each vial is crafted using premium-quality, chemically-resistant glass that ensures excellent durability and protection, and providing unparalleled clarity.

Range of Sizes

We understand that every project is unique. That's why our new range is available in a variety of sizes, from 2 ml to 100 ml, catering to all your specific storage needs.

Environmental Responsibility

We're committed to protecting the planet. Our unique and innovative packaging is 100% recyclable and eco-friendly, without compromising on quality or durability.

If you require any more information on this product, please contact our customer service team.