

PRODUCT DATA SHEET



CV0024

Rolled Rim Vials, With Unattached Clip On / Snap On Closure



Product comes with unattached PE clip on / snap on closure

Specifications

Vial Material	Glass Borosilicate 3 (Autoclavable)
Capacity	14ml
Outer Diameter	25.5mm
Height	48mm
Closure Material	Polyethylene (PE)
Closure Type	Clip On / Snap
Closure Liner	N/A
Pack Quantity	100
Packaging	Cardboard Box & Plastic Bag
Pack Length	420mm
Pack Width	240mm
Pack Height	54mm
Pack Weight	1.8Kg
Case Quantity	10packs (1000vials)
Case Length	500mm
Case Width	440mm
Case Height	290mm
Case Weight	18.2Kg
Pallet Quantity	200packs (20000Vials)

What's in the box?



Streamline your workflow and simplify your storage with our pack of 100 vials and caps. They neatly stack and nest together, maximizing your space and allowing convenient access whenever needed.

Unrivalled Quality

Each vial is crafted using premium-quality, chemically-resistant glass that ensures excellent durability and protection, and providing unparalleled clarity.

Range of Sizes

We understand that every project is unique. That's why our new range is available in a variety of sizes, from 7ml to 30ml, catering to all your specific storage needs.

Environmental Responsibility

We're committed to protecting the planet. Our unique and innovative packaging is 100% recyclable and eco-friendly, without compromising on quality or durability.

If you require any more information on this product, please contact our customer service team.