# **PRODUCT DATA SHEET**



#### <u>CV0046</u>

Universal Vial, Tubular, White Polypropylene Closure



Product comes with unattached polypropylene closure

## **Specifications**

Vial Material	Glass, Soda (Not recommended for autoclave)
Capacity	28ml
Outer Diameter	27.5mm
Height	82.5mm
Closure Type	Screw, White Polypropylene
Closure Liner	Unlined
Pack Quantity	100
Packaging	Cardboard Box
Pack Length	448mm
Pack Width	260mm
Pack Height	870mm
Pack Weight	2.6Kg
Case Quantity	5packs (500vials)
Case Length	470mm
Case Width	280mm
Case Height	460mm
Case Weight	13.8Kg
Pallet Quantity	160packs (16000Vials)



The Bijou Vial, tubular vial designed for the secure storage and transport of liquids and small samples. Each vial is equipped with a polypropylene closure, maintaining the integrity and sterility of the contents.

### Available Sizes

- **7ml Vial CV0045:** Ideal for smaller sample volumes, the 7ml Bijou Vial is compact and convenient for a variety of laboratory applications.
- **28ml Vial CV0046:** The 28ml Universal Vial offers a larger capacity, suitable for storing larger quantities of samples or reagents.

### <u>Key Features</u>

- **Tubular Construction:** The vials are produced to precise specifications, ensuring uniformity and high performance in demanding laboratory environments.
- **Polypropylene Closure:** Each vial has and unattached polypropylene closure, providing a secure seal that protects the contents from contamination and evaporation.
- **Versatility:** Suitable for a wide range of applications including clinical, pharmaceutical, and research laboratories.

### **Environmental Responsibility**

We're committed to protecting the planet. Our unique and innovative packaging is 100% recyclable and eco-friendly, without compromising on quality or durability.

If you require any more information on this product, please contact our customer service team.