

Contour - 1000 Media Bottle

-  1000ml Capacity, moulded 360° graduations
-  Diameter 96.5mm +/- 1.5
-  Overall height 241mm +/- 1.5
-  Amber Type III Glass
-  Autoclavable
-  Recyclable



An opportunity to clean up!

The Price is Favourable, and the Environmental Impact Substantial.



Contour - 1000 Media Bottle 

The price is favourable and the environmental impact substantial



- Avoid Lab Specification overkill
- Leakproof GL45 cap to fit media bottle attachments
- UV Protection
- Profile 360° graduations serve as extra grip
- Made with over 70% Recycled Material
- USP 660 Compliant



Also known as 'The Bursar's Favourite', this highly economical, UV protected media bottle alternative will enhance your Sustainability status, as well as help you comply with budget tightening.

The ergonomic 360° graduations offer enhanced grip to reduce accidents

- Easy disposal in any bottle bank improves sustainability and recycling (the vessel itself will last for decades, there is no need to throw it away)
- High recycled material content is an environmental plus
- The GL45 cap is fully compatible with other media bottles
- Answers your budget requirements.

Autoclave approved: re-use over and over again

There are occasions when you require the very best in terms of Thermal Shock and Hydrolytic Resistance. However, where a compromise is possible, and it often is, then the Contour-1000 will offer great savings and only relatively small differences in performance. TO AVOID LAB CONFUSION WITH BOROSILICATE 3.3 BOTTLES, THE CONTOUR-1000 HAS PROFILED GRADUATIONS TO CLEARLY DIFFERENTIATE IT AS A TYPE III USP 660 MATERIAL.

This bottle is autoclavable, so can be sterilised and used over and over again; we have performed 50 back to back autoclave cycles. It is recyclable, and manufactured from over 70% recycled material. It is also very well priced.

Our point is that judicious lab management avoids overkill, and there are a plethora of instances where the Contour-1000 will easily suffice, leaving you free to book the string quintet as opposed to the whole orchestra. If you are in any doubt, just ask the Bursar!



No Room for Overkill